


## TekBIO ELECTROPHORESIS DC power supplies

TekBIO electrophoresis Dc power supplies are long lasting and built using quality components, It can supply pure DC resulting in sharper gel bands

### TekBIO Tapped D.C Voltage power supply-50/100V


Tapped D.C voltage power supply 50 and 100V @100mA. 1 set 4 mm banana sockets and user manual



Pro No	FIXI-50/100
Body construction	MS powder coated
Input Voltage	230 V AC +/- 10%, 50 Hz
No of Outputs	1 No
Output voltages	D.C- 50 and 100V
Output current	100mA
Voltage Display	LED bulb

### TekBIO FIXI- 100 Power Supply

Tapped-Fixed D.C voltage power supply 50, 100, 150 and 200V @ 100mA. Output-1 No and user manual



Pro No	FIXI-100
Body construction	MS powder coated
Input Voltage	230 V AC +/- 10%, 50 Hz
No of Outputs	1 No
Output voltages	D.C- 50, 100, 150 & 200V
Output current	100mA
Voltage Display	LED bulb

## TekBIO FIXI-200 Power Supply

Tapped - Fixed D.C voltage power supply 50, 100, 150 and 200V @ 200mA.  
Output - 1 No and user manual



Pro No	FIXI-200
Body construction	MS powder coated
Input Voltage	230 V AC +/- 10%, 50 Hz
No of Outputs	1 No
Output voltages	D.C- 50, 100, 150 & 200V
Output current	200mA
Voltage Display	LED bulb

## TekBIO DIGI-100 Digital Power supply with Digital Timer-TModel

Digital D.C Power supply, 0-300V, 0-100mA, Constant Voltage, Constant Current, With Digital display and Digital Timer, output 1-No and user manual



Pro No	DIGI-100T
Body construction	MS powder coated
Input Voltage	230 V AC +/- 10%, 50 Hz
No of Outputs	2 No
Output voltages	0-300V insteps of 1 V
Output current	0-100mA insteps of 1 mA
Voltage Display	Digital LED display
Digital Timer	YES

## TekBIO DIGI-200 Digital Power supply with Digital Timer-TModel


Digital D.C Power supply, 0-300V, 0-200mA, Constant Voltage, Constant Current, With Digital display and Digital Timer, output 1-No and user manual



Pro No	DIGI-200T
Body construction	MS powder coated
Input Voltage	230 V AC +/- 10%, 50 Hz
No of Outputs	2 No
Output voltages	0-300V insteps of 1 V
Output current	0-200mA insteps of 1 mA
Voltage Display	Digital LED display
Digital Timer	YES

## TekBIO DIGI-500 Digital Power supply with Digital Timer-SMPS-Model

Digital D.C Power supply, 0-500V, 0-500mA, Constant Voltage, Constant Current, With Digital display and Digital Timer, output - 4 No and user manual



Pro No	DIGI-500T
Body construction	MS powder coated
Input Voltage	230 V AC +/- 10%, 50 Hz
No of Outputs	4 No
Output voltages	0-500V insteps of 1 V
Output current	0-500mA insteps of 1 mA
Voltage Display	Digital LED display
Digital Timer	YES