

TekBIO Vortex Mixer



A vortex mixer, or vortexer, is a simple device used commonly in laboratories to mix small vials of liquid. When a test tube or other appropriate container is pressed into the rubber cup (or touched to its edge) the motion is transmitted to the liquid inside and a vortex is created.

Product Details:

Variable Speed Setting	2500 rpm
Operating Modes	Continuous
Sample Capacity/Format	48 wells
Volume Thermal Block Sample	0.1 ml and 0.2 ml

*Product images are for illustrative purposes only. Actual product may vary from the Image.

